




RADIX

New RADIX systems are intended for RO or Grade 2 water supply central laboratory house water distribution systems. It can be connected to atmospheric tanks with level sensors.

DESCRIPTION

	RADIX 60 RO	RADIX 100 RO	RADIX 30 EDI	RADIX 60 EDI
Model	RA1001	RA1003	RA2101	RA2103
Capacity, l/H	60	100	30	60
Average ion rejection rate	99%	99%	-	-
Output water conductivity	-	-	0.2 μ S	0.2 μ S
Display	Controller with color display			
Dimensions (WxDxH), cm	33x65x72	33x65x72	33x65x72	33x65x72
Weight, kg	42	56	50	56

FLOW DIAGRAM



Dzerbenes 27, Riga
LATVIA, LV-1006
Tel.: +371 67551894
Fax.: +371 67551976
info@adrona.lv
www.adrona.lv